Multisil-Epithetik soft-form

Multisil-Epithetik hard-form

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Multisil-Epithetik soft-form

Multisil-Epithetik hard-form

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Indication range

Multisil-Epithetik soft-form and Multisil-Epithetik hard-form on the basis of addition-linking vinyl poly- dimethylsiloxanes are used for oral and extraoral impressions in the field of epithetics. Multisil-Epithetik soft-form is particularly suitable to reproduce undercut areas and problematic zones.

Multisil-Epithetik hard-form is used to prepare flat impressions and as reinforcing material in combination with Multisil-Epithetik soft-form.

2. Indication

Oral and extraoral impressions in the field of epithetics.

3. Safety recommendations and hazard warnings

After skin contact, wash with copious amount of water and soap. After eye contact, rinse with copious amount of water. Please request the safety datasheets of the product separately.

3.1 Symbols

	Can be used until YYYY-MM
	Lot number

4. Storage and durability information

4.1 Storage

Multisil Epithetik soft-form and hard-form must be stored under dry and cool conditions (not above 28° C). Use the original cap to seal the cartridge after use; protect against sunlight.

4.2 Durability

In the original, sealed condition - see label.

5. Technical data

Processing time: 1 to 11/2 minutes at 23 °C, at higher temperature the processing time is reduced, at lower temperatures extended.

Can be removed from the mould after 7 to 10 minutes at 23 °C (after approx. 5 minutes at 37° C).

Final hardness of Multisil-Epithetik soft-form: 25 Shore A

Final hardness of Multisil-Epithetik hard-form: 45 Shore A

6. Interaction with other materials

The highly reactive catalyst that is contained in addition-linking silicones reacts highly sensitive to certain substances („catalyst poisons“). The most important catalyst poisons are organic sulphur compounds, amines and amides and compounds containing heavy metals such as lead, tin, cadmium or mercury. Contact of Multisil-Soft which has not hardened with such compounds (e.g. with latex gloves, partly also with condensation-linking silicones) results in damage to the polymerization behaviour (inhibition, retarded hardening, formation of a grease layer) and must be avoided.

7. Processing

Insert the Multisil-Epithetik soft-form or hard-form double cartridge which is not compressed (pay attention to the notch). Turn the protective cap to the left to screw it off the cartridge. Place on the mixing cannula (pay attention to the notch) and fix by turning to the right. By exerting uniform pressure onto the manual lever the material is fed into the mixing cannula and forcefully mixed. If new cartridges are used, discard a strip of a length of 1 cm, since due to manufacturing-technical reasons there may be slight differences in the quantity of material filled into the cartridge halves.

Multisil-Epithetik soft-form or hard-form are applied uniformly and without the formation of bubbles. Bubble-free application can be achieved if the area to be reproduced is coated by exerting uniform pressure and without stopping the application process.

Once the material has been applied, remove the mixing cannula and place the original cap on again. The safety cap avoids contamination of the discharge holes of the cartridge.

Due to its low Shore A hardness, Multisil-Epithetik soft-form needs to be coated with Multisil-Epithetik hard-form; large surfaces must be stabilized using small wooden sticks.

8. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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